



Product Name: Recombinant Human POLR1G Protein
Catalog #: PEH2720



Summary

Name	RPA34_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15446
Host	E.coli
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEEPQAGDAA RFSCPPNFTA KPPASESPRF SLEALTGPDT ELWLIQAPAD FAPECFNGRH VPLSGSQIVK GKLAGKRHRY RVLSSCPQAG EATLLAPSTE AGGGLTCASA PQGTLRILEG PQQSLSGSPL QPIPASPPPQ IPPGLRPRFC AFGGNPPVTG PRSALAPNLL TSGKKKKEMQ VTEAPVTQEA VNGHGALEVD MALGSPEMDV RKKKKKKKNQQ LKEPEAAGPV GTEPTVETLE PLGVLPSTT KKRKKPKGKE TFEPEDKTVK QEQUINTEPLE DTVLSPTKKR KRQKGTEGME PEEGVTVESQ PQVKVEPLEE AIPLPPTKKR KKEKGQMAMM EPGTEAMEPV EPEMKPLESP GGTMAPQQPE GAKPQAQAAL AAPKKKTKKE KQQDATVEPE TEVVGPELPD DLEPQAPTS TTKKKKKKKER GHTVTEPIQP LEPELPGEGQ PEARATPGST KKRKKQSQES RMPETVPQEE MPGPPLNSES GEEAPTGRDK KRKQQQQQPV

image



Background

Product Name: Recombinant Human POLR1G Protein
Catalog #: PEH2720



Alternative Names

DNA-directed RNA polymerase I subunit RPA34 (A34.5) (Antisense to ERCC-1 protein) (ASE-1) (CD3-epsilon-associated protein) (CD3E-associated protein) (DNA-directed RNA polymerase I subunit G) (RNA polymerase I-associated factor PAF49)

Note

For Research Use Only , Not for Diagnostic Use