


Product Name: Recombinant Human FLRT2 Protein
Catalog #: PEH2722



Summary

Name	FLRT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43155
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RCDRNFVYCN ERSLTSVPLG IPEGVTVLVL HNNQINNAGF PAELHNVQSV HTVYLYGNQL DEFPMNLPKN VRVLHLQENN IQTISRAALA QLLKLEELHL DDNSISTVGV EDGAFREAIS LKLLFLSKNH LSSVPVGLPV DLQELRVDEN RIAVIDMAF QNLTSLERLI VDGNLLTNKG IAEGTFSHLT

image



Background

Alternative Names	Leucine-rich repeat transmembrane protein FLRT2 (Fibronectin-like domain-containing leucine-rich transmembrane protein 2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FLRT2 Protein
Catalog #: PEH2722

