


**Product Name: Recombinant Human FLRT2 Protein**  
**Catalog #: PEH2723**



## Summary

|                              |                                                                                                                                                                                                                          |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | FLRT2_HUMAN                                                                                                                                                                                                              |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                      |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                         |
| <b>Tag</b>                   | N-SUMO & C-His                                                                                                                                                                                                           |
| <b>Accession #</b>           | O43155                                                                                                                                                                                                                   |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                   |
| <b>Species</b>               | Human                                                                                                                                                                                                                    |
| <b>Molecular Mass</b>        | 35 kDa                                                                                                                                                                                                                   |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                       |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.           |
| <b>Sequence</b>              | RCDRNFVYCN ERSLTSVPLG IPEGVTVLVL HNNQINNAGF PAELHNVQSV<br>HTVYLYGNQL DEFPMNLPKN VRVLHLQENN IQTISRAALA QLLKLEELHL<br>DDNSISTVGV EDGAFREAIS LKLLFLSKNH LSSVPVGLPV DLQELRVDEN<br>RIAVIDMAF QNLTSLERLI VDGNLLTNKG IAEGTFSHLT |

## image



## Background

|                          |                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Leucine-rich repeat transmembrane protein FLRT2 (Fibronectin-like domain-containing leucine-rich transmembrane protein 2) |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------|

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human FLRT2 Protein**  
**Catalog #: PEH2723**

---

