

Product Name: Recombinant Human NEMF Protein
Catalog #: PEH2730



Summary

Name	NEMF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60524
Host	E.coli
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DIEKVLVSLQ KAEDYMKTTS NFSGKGYYIQ KREIKPSLEA DKPVEDILTY EEFHPFLFSQ HSQCPYIEFE SFDKAVDEFY SKIEGQKIDL KALQQEKQAL KKLDNVRKDH ENRLEALQQA QEIDKLKGEL IEMNLQIVDR AIQVVRSA NQIDWTEIGL IVKEAQAQGD PVASAIKELK LQTNHVTMLL RNPYLLSEEE DDDVDGDVNV EKNETEPPKG KKKKQKNKQL QKPQKNKPLL VDVDLSLSAY ANAKKYYDHK RYAAKKTQKT VEAAEKAFKS AEKKTQTLK EVQTVTSIQK ARKVYWFEKF LWFISSENYL IIGGRDQQQN EIIVKRYLTP GDIYVHADLH GATSCVIKNP TGEPIPPRTL TEAGTMALCY SAAWDARVIT SAWWVYHHQV SKTAPTGEYL TTGSFMIRGK KNFLPPSYLM MGFSFLFKVD ESCVWRHQGE RKVRVQDEDM ETLASCTSEL ISEEMEQLDG GDTSSDEDKE EHETPVEVEL MTQVDQEDIT LQSGRDELNE ELIQEESSED EGEYEEVRKD QDSVGEMKDE GEETLNYPDT TIDLSHLQPQ

image



Product Name: Recombinant Human NEMF Protein
Catalog #: PEH2730



Background

Alternative Names

Nuclear export mediator factor NEMF (Antigen NY-CO-1) (Serologically defined colon cancer antigen 1)

Note

For Research Use Only , Not for Diagnostic Use