


Product Name: Recombinant Human JAK2 Protein
Catalog #: PEH2731



Summary

Name	JAK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60674
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PTQFEERHLK FLQQLGKGNF GSVEMCRYDP LQDNTGEVVA VKKLQHSTEE HLRDFEREIE ILKSLQHDNI VKYKGVCYSA GRRNLKLIME YLPYGSLRDY LQKHKERIDH IKLLQYTSQI CKGMEYLGTK RYIHRDLATR NILVENENRV KIGDFGLTKV LPQDKEYYKV KEPGESPIFW YAPESLTESK FSVASDVWSF GVVLYELFTY IEKSKSPPAE FMRMIGNDKQ GQMIVFHLIE LLKNNGRLPR PDGCPDEIYM IMTECWNNNV NQRPSFRDLA LRVDQIRDNM

image



Background

Alternative Names	Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2)
--------------------------	---

Product Name: Recombinant Human JAK2 Protein
Catalog #: PEH2731



Note

For Research Use Only , Not for Diagnostic Use