


Product Name: Recombinant Human WDR1 Protein
Catalog #: PEH2737



Summary

Name	WDR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	O75083
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLPQVERGVS KIIGDPPKGN NFLYTNGKCV ILRNIDNPAL ADIYTEHAHQ VVVAKYAPSG FYIASGDVSG KLRIWDTTQK EHLLKYEQP FAGKIKDIAW TEDSKRIAVV GEGREKFGAV FLWDSGSSVG EITGHNKVIN SVDIKQSRPY RLATGSDDNC AAFEGPPFK FKFTIGDHSR FVNCVRFSPD GNRFATASAD GQIYIDGKT GEKVCALGGS KAHDGGIYAI SWSPDSTHLL SASGDKTSKI

image



Background

Alternative Names WD repeat-containing protein 1 (Actin-interacting protein 1) (AIP1) (NORI-1)

Note

Product Name: Recombinant Human WDR1 Protein
Catalog #: PEH2737



For Research Use Only , Not for Diagnostic Use