


Product Name: Recombinant Human WDR1 Protein
Catalog #: PEH2738



Summary

Name	WDR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75083
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLRIWDTTQK EHLKYEQP FAGIKDIAW TEDSKRIAVV GEGREKFGAV FLWDSGSSVG EITGHNVIN SVDIKQSRPY RLATGSDDNC AAFEGPPFK FKFTIGDHSR FVNCVRFSPD GNRFATASAD GQIYIDGKT GEKVCALGGS KAHDGGIYAI SWSPDSTHLL SASGDKTSKI WDVSVNSVVS TFPMGSTVLD QQLGCLWQKD HLLSVLSGY INYLDNRNPS KPLHVIKGHS KSIQCLTVHK NGGKSYIYSG SHDGHINYWD SETGENDSFA GKGHTNQVSR

image



Background

Alternative Names WD repeat-containing protein 1 (Actin-interacting protein 1) (AIP1) (NORI-1)

Product Name: Recombinant Human WDR1 Protein
Catalog #: PEH2738



Note

For Research Use Only , Not for Diagnostic Use