



Product Name: Recombinant Human CILP Protein
Catalog #: PEH2741



Summary

Name	CILP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75339
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NISTATAAQT DLNFINDEGD TFPLRTYGMF SVDFRDEVTS EPLNAGKVKV HLDSTQVKMP EHISTVKLWS LNPDTGLWEE

image



Background

Alternative Names	Cartilage intermediate layer protein 1 (CILP-1) (Cartilage intermediate-layer protein) [Cleaved into: Cartilage intermediate layer protein 1 C1; Cartilage intermediate layer protein 1 C2]
--------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use