



Product Name: Recombinant Human DKK1 Protein
Catalog #: PEH2750



Summary

Name	DKK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O94907
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TLNSVLNSNA IKNLPPPLGG AAGHPGSAVS AAPGILYPGG NKYQTIDNYQ PYPCAEEDEC GTDEYCASPT RGGDAGVQIC LACRKRKRRC MRHAMCCPGN YCKNGICVSS DQNHFRGEIE ETITESFGND HSTLDGYSRR TTLSSKMYHT KGQEGSVCLR SSDCASGLCC ARHFWSKICK PVLKEGQVCT KHRRKGS HGL EIFQRCYCGE GLSCRIQKDH HQASNSSLRH TCQRH

image



Background

Alternative Names Dickkopf-related protein 1 (Dickkopf-1) (Dkk-1) (hDkk-1) (SK)

Note

Product Name: Recombinant Human DKK1 Protein
Catalog #: PEH2750



For Research Use Only , Not for Diagnostic Use