


Product Name: Recombinant Human ABCA1 Protein
Catalog #: PEH2752



Summary

Name	ABCA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95477
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HESDTLTIDV SAISNLIRKH VSEARLVEDI GHETYVLPY EAAKEGAFVE LFHEIDDRLS DLGISSYGIS ETTLEEFLK VAEESGVDAE TSDGTLPARR NRRAFGDKQS

image



Background

Alternative Names	Phospholipid-transporting ATPase ABCA1 (EC 7.6.2.1) (ATP-binding cassette sub-family A member 1) (ATP-binding cassette transporter 1) (ABC-1) (ATP-binding cassette 1) (Cholesterol efflux regulatory protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ABCA1 Protein
Catalog #: PEH2752

