


Product Name: Recombinant Human F13A1 Protein
Catalog #: PEH2757



Summary

Name	F13A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P00488
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FVFAEVNSDL IYITAKKDGT HVVENVDATH IGKLIVTKQI GGDGMMMDITD TYKFQEGQEE ERLALETALM YGAKKPLNTE GVMKRSRNVD

image



Background

Alternative Names	Coagulation factor XIII A chain (Coagulation factor XIIIa) (EC 2.3.2.13) (Protein-glutamine gamma-glutamyltransferase A chain) (Transglutaminase A chain)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use