

Product Name: Recombinant Human EGF Protein
Catalog #: PEH2771



Summary

Name	EGF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P01133
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLLSKNPKNP YEESSRDVRS RRPADTEDGM SSCPQPWFVV IKEHQDLKNG GQPVAGEDGQ AADGSMQPTS WRQEPQLCGM

image



Background

Alternative Names	Pro-epidermal growth factor (EGF) [Cleaved into: Epidermal growth factor (Urogastrone)]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use