


Product Name: Recombinant Human MBP Protein
Catalog #: PEH2792



Summary

Name	MBP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P02686
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GREDNTFKDR PSEDELQTI QEDSAATSES LDVMASQKRP SQRHGSKYLA TASTMDHARH GFLPRHRDTG ILDSIGRFFG GDRGAPKRGs

image



Background

Alternative Names	Myelin basic protein (MBP) (Myelin A1 protein) (Myelin membrane encephalitogenic protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use