


Product Name: Recombinant Human MMP1 Protein
Catalog #: PEH2804



Summary

Name	MMP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P03956
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QTHLTYRIEN YTPDLPRADV DHAIEKAFQL WSNVTPLTFT KVSEGGQADIM ISFVRGDHRD NSPFDGPGGN LAHAFQPGPG IGGDAHFDDED ERWTNNFREY NLHRVAAHEL GHSLGLSHST DIGALMYPSTY TFGSDVQLAQ DDIDGIQAIY GRSQNPVQPI GPQTPKACDS KLTFDAITTI

image



Background

Alternative Names	Interstitial collagenase (EC 3.4.24.7) (Fibroblast collagenase) (Matrix metalloproteinase-1) (MMP-1) [Cleaved into: 22 kDa interstitial collagenase; 27 kDa interstitial collagenase]
--------------------------	---

Note

Product Name: Recombinant Human MMP1 Protein
Catalog #: PEH2804



For Research Use Only , Not for Diagnostic Use