


Product Name: Recombinant Human PROC Protein
Catalog #: PEH2807



Summary

Name	PROC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P04070
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RHSSLERECI EEICDFEAK EIFQNVDDL AFWSKHVDGD QCLVPLEHP CASLCCGHGT CIDGIGSFSC DCRSGWEGRF CQREVSFLNC SLDNNGCTHY CLEEVGWRRRC SCAPGYKLGD DLLQCHPAVK FPCGRPWKRM

image



Background

Alternative Names Vitamin K-dependent protein C (EC 3.4.21.69) (Anticoagulant protein C) (Autoprothrombin IIA) (Blood coagulation factor XIV) [Cleaved into: Vitamin K-dependent protein C light chain; Vitamin K-dependent protein C heavy chain; Activation peptide]

Note

Product Name: Recombinant Human PROC Protein
Catalog #: PEH2807



For Research Use Only , Not for Diagnostic Use