


Product Name: Recombinant Human VWF Protein
Catalog #: PEH2821



Summary

Name	VWF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P04275
Host	E.coli
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRNSMVL DVA FVLEGS D KIG EADFNR SKEF MEEVIQR MDV GQDSIHVT VL QYSYMVT VEY PFSEAQ SKGD ILQRVRE IRY QGGNRT NTGL ALRYLSDH SF LVSQGDRE QA PNLVYMVT GN PASDEIK RLP GDIQVVP IGV GPNANVQE LE RIGWPNAP IL IQDFETL PRE APDLVLQR CC SG EGLQIPT L SPAPDCSQ PL DVILLLDG SS SFPASYF DEM KSF AKAFISK ANIGPRLT QV SVLQYGSITT IDVPWNV VPE KAHLLSLVD V MQREGGPSQ I GDALGF AVRY LTSEM HGARP GASKAVVIL V TDVSVDSV DA AADAARSNR V TVFPIGIG DR YDAAQLRILA GPAGDSNVVK LQRIEDLPTM VTLGNSFLHK LCSGFV RICH DEDGNEKRPG DVWTLPDQCH

image



Background

Product Name: Recombinant Human VWF Protein
Catalog #: PEH2821



Alternative Names

von Willebrand factor (vWF) [Cleaved into: von Willebrand antigen 2 (von Willebrand antigen II)]

Note

For Research Use Only , Not for Diagnostic Use