


Product Name: Recombinant Human PRM2 Protein
Catalog #: PEH2823



Summary

Name	PRM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P04554
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERSHEVYRQQ LHGQEQGHHG QEEQGLSPEH VEVYERTHGQ SHYRRRHCSR RRLHRIHRRQ HRSCRRRKRR

image



Background

Alternative Names	Protamine-2 (Sperm histone P2) (Sperm protamine P2) [Cleaved into: Basic nuclear protein HPI1; Basic nuclear protein HPI2; Basic nuclear protein HPS1; Basic nuclear protein HPS2; Sperm histone HP4 (Sperm protamine P4); Sperm histone HP2 (Sperm protamine P2) (P2'); Sperm histone HP3 (P2'') (Sperm protamine P3)]
--------------------------	---

Note

Product Name: Recombinant Human PRM2 Protein
Catalog #: PEH2823



For Research Use Only , Not for Diagnostic Use