



Product Name: Recombinant Human CYP1A1 Protein
Catalog #: PEH2827



Summary

Name	CP1A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P04798
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HQELLSLVNL NNNFGEVVGSGNPADFIPILRYLPNPSLNAFKDLNEKFYS FMQKMOVKEHYKTFEKGHIRDITDSLIEHCQEKQLDENANV

image



Background

Alternative Names	Cytochrome P450 1A1 (CYP1A1) (EC 1.14.14.1) (Cytochrome P450 form 6) (Cytochrome P450-C) (Cytochrome P450-P1) (Hydroperoxy icosatetraenoate dehydratase) (EC 4.2.1.152)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use