


**Product Name: Recombinant Human ALPP Protein**  
**Catalog #: PEH2835**



## Summary

<b>Name</b>	PPB1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P05187
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	53 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	PDFWNREAAE ALGAAKKLQP AQTAAKNLII FLGDGMGVST VTAARILKGGQ KKDKLGPEIP LAMDRFPYVA LSKTYNVDKH VPDSGATATA YLCGVKGNFQ TIGLSAAARF NQCNTTRGNE VISVMNRAKK AGKSVGVT TTRVQHASPAG TYAHTVNRNW YSDADVPASA RQEGCQDIAT QLISNMDIDV ILGGGRKYM RMGTPDPEYP DDYSQGGTRL DGKNLVQEWL AKRQGARYVW NRTELMQASL DPSVTHLMGL FEPGDMKYEI HRDSTLDPSL MEMTEAALRL LSRNPRGFFL FVEGGRIDHG HHESRAYRAL TETIMFDDAI ERAGQLTSEE DTLSLVTADH SHVFSFGGYP LRGSSIFGLA PGKARDRKAY TVLLYGNGPG YVLKDGARPD VTESEGSPE YRQSAVPLD EETHAGEDVA VFARGPQ AHL VHGVQEQTFI AHVMAFAACL EPYTACDLAP

## image



## Background

**Product Name: Recombinant Human ALPP Protein**  
**Catalog #: PEH2835**



---

**Alternative Names**

Alkaline phosphatase, placental type (EC 3.1.3.1) (Alkaline phosphatase Regan isozyme) (Placental alkaline phosphatase 1) (PLAP-1)

**Note**

For Research Use Only , Not for Diagnostic Use