


Product Name: Recombinant Human KRT8 Protein
Catalog #: PEH2840



Summary

Name	K2C8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P05787
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KFASFIDKVR FLEQQNKMLE TKWSLLQQQK TARSNMDNMF ESYINNLRRQ LETLGQEKLK LEAELGNMQG LVEDFKNKYE DEINKRTEME NEFVLIKKDV DEAYMNKVEL ESRLEGLTDE INFLRQLYEE EIRELQSQIS DTSVVLMDN SRSLDMDSII AEVKAQYEDI ANRSRAEAE MYQIKYEELQ SLAGKHGDDL RRTKTEISEM NRNISRQAIE IGLKQQRAS LEAAIADAEQ RGELAIKDAN AKLSELEAAL QRAKQDMARQ LREYQELMNV KLALDIEIAT YRKLLEGEES RLESGMQNMS IHTKTTSGYA GGLSSAYGGL TSPGLSYSLG SSFGSGAGSS SFRSTSSRA VVVKKIETRD GKLVSSESDV

image



Background

Alternative Names Keratin, type II cytoskeletal 8 (Cytokeratin-8) (CK-8) (Keratin-8) (K8) (Type-II)

Product Name: Recombinant Human KRT8 Protein
Catalog #: PEH2840



keratin Kb8)

Note

For Research Use Only , Not for Diagnostic Use