


Product Name: Recombinant Human INSR Protein
Catalog #: PEH2841



Summary

Name	INSR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P06213
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DERTYGAKS DIIVVQTDAT NPSVPLDPIS VSNSSSQIIL KWKPPSDPNG NITHYLVFWE RQAEDSELF LDYCLKGLKL PSRTWSPPF SEDSQKHNQS EYEDSAGECC SCPKTD SQIL KELEESSFRK TFEDYLHNVV FVPRKTSSGT

image



Background

Alternative Names	Insulin receptor (IR) (EC 2.7.10.1) (CD antigen CD220) [Cleaved into: Insulin receptor subunit alpha; Insulin receptor subunit beta]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use