


Product Name: Recombinant Human PTMA Protein
Catalog #: PEH2845



Summary

Name	PTMA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P06454
Host	E.coli
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSDAAVDTSS EITTKDLKEK KEVVEEAENG RDAPANGNAE NEENGEQEAD NEVDEEEEEEG GEEEEEEEEEG DGEEEDGDED EEAESATGKR AAEDDEDDDV DTKKQKTDED D

image



Background

Alternative Names	Prothymosin alpha [Cleaved into: Prothymosin alpha, N-terminally processed; Thymosin alpha-1]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use