


**Product Name: Recombinant Human NEFM Protein**  
**Catalog #: PEH2854**



## Summary

<b>Name</b>	NFM_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P07197
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	66 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	EKEQLQGLND RFAGYIEKVH YLEQQNKEIE AEIQALRQKQ ASHAQLGDAY DQEIRELRAT LEMVNHEKAQ VQLDSHLEE DIHRLKERFE EEARLRDDTE AAIRALRKDI EEASLVKVEL DKKVQSLQDE VAFLRSNHEE EVADLLAQIQ ASHITVERKD YLKTDISTAL KEIRSQLESH SDQNMHQAEW WFKCRYAKLT EAAEQNKEAI RSAKEEIAEY RRQLQSKSIE LESVRGTKES LERQLSDIEE RHNHDLSSYQ DTIQQLENEL RGTKWEMARH LREYQDLLNV KMALDIEIAA YRKLEGEET RFSTFAGSIT GLYTHRPPI TISSKIQKPK VEAPKLVQH KFVEEIIIEET KVEDEKSEME EALTAITEEL AVSMKEEKKE AAEEKEEPE AEEEEVAACK SPVKATAPEV KEEGEKEEEE EGQEEEEED EGAKSDQAEW GGSEKEGSSE KEEGEQEEGE TEAAEGEEA EAKEEKKVEE KSEEVATKEE LVADAKVEKP EKAKSPVPS PVEEKGKSPV PKSPVEEKGK SPVPSVVEE

## image



## Background

**Product Name: Recombinant Human NEFM Protein**  
**Catalog #: PEH2854**



---

**Alternative Names**

Neurofilament medium polypeptide (NF-M) (160 kDa neurofilament protein)  
(Neurofilament 3) (Neurofilament triplet M protein)

**Note**

For Research Use Only , Not for Diagnostic Use