


Summary

Name	NFM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P07197
Host	E.coli
Species	Human
Molecular Mass	62 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKEQLQGLND RFAGYIEKVH YLEQQNKEIE AEIQALRQKQ ASHAQLGDAY DQEIRELRAT LEMVNHEKAQ VQLDSDHLEE DIHRLKERFE EEARLRDDTE AAIRALRKDI EEASLVKVEL DKKVQSLQDE VAFLRSNHEE EVADLLAQIQ ASHITVERKD YLKTDISTAL KEIRSQLESH SDQNMHQAEW WFKCRYAKLT EAAEQNKEAI RSAKEEIAEY RRQLQSKSIE LESVRGTKES LERQLSDIEE RHNHDLSSYQ DTIQQLENEL RGTKWEMARH LREYQDLLNV KMALDIEIAA YRKLEGEET RFSTFAGSIT GPLYTHRPPi TISSKIQKPK VEAPKLKVQH KFVEEIIIEE KVEDEKSEME EALTAITEEL AVSMKEEKKE AAEEKEEPE AEEEEVAACK SPVKATAPEV KEEGEKEEEE EGQEEEEED EGAKSDQAEW GGSEKEGSSE KEEGEQEEGE TEAEAEGEEA EAKEEKKVEE KSEEVATKEE LVADAKVEKP EKAKSPVPKS

image



Background

Product Name: Recombinant Human NEFM Protein
Catalog #: PEH2855



Alternative Names

Neurofilament medium polypeptide (NF-M) (160 kDa neurofilament protein)
(Neurofilament 3) (Neurofilament triplet M protein)

Note

For Research Use Only , Not for Diagnostic Use