


**Product Name: Recombinant Human TUBB Protein**  
**Catalog #: PEH2867**



## Summary

|                              |   |
|------------------------------|---|
| <b>Name</b>                  | TBB5_HUMAN  |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE   |
| <b>Endotoxin level</b>       | Please contact with the lab for this information  |
| <b>Tag</b>                   | C-His   |
| <b>Accession #</b>           | P07437  |
| <b>Host</b>                  | E.coli  |
| <b>Species</b>               | Human   |
| <b>Molecular Mass</b>        | 50 kDa  |
| <b>Buffer</b>                | Phosphate buffered saline   |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.  |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.  |
| <b>Sequence</b>              | QCGNQIGAKF WEVISDEHGI DPTGTYHGDS DLQLDRISVY YNEATGGKYV<br>PRAILVDLEP GTMDSVRS GP FGQIFRPDNF VFGQSGAGNN WAKGHYTEGA<br>ELVDSVLDVV RKEAESCDC L QGFQLTHSLG GGTGSGMGTL LISKIR E EYP<br>DRIMNTFSVV PSPKVS DTVV EPYNATLSVH QLVENTDETY CIDNEALYDI<br>CFRTLKLTTP TYGDLNHLVS ATMSGVTTCL RFPGQLNADL RKLAVNMV PF<br>PRLHFFMPGF APLTSRGSQQ YRALTVP ELT QQVFD AKNMM AACDP RHGRY<br>LTVA AVFRGR MSMKEVDEQM LNVQNKSSY FVEWIPNNVK TAVCDIPPRG<br>LKMAVTFIGN STAIQELFKR ISEQTAMFR RKAFLHWYTG EGMD EMEFTE<br>AESNMNDLVS EYQQYQDATA EEEEDFGEEA |

## image



## Background

**Product Name: Recombinant Human TUBB Protein**  
**Catalog #: PEH2867**



---

**Alternative Names**

Tubulin beta chain (Tubulin beta-5 chain)

**Note**

For Research Use Only , Not for Diagnostic Use