


Product Name: Recombinant Human DCN Protein
Catalog #: PEH2868



Summary

Name	PGS2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P07585
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPFQQRGLFD FMLEDEASGI GPEVPDDRDF EPSLGPVCPF RCQCHLRVVQ CSDLGLDKVP KDLPPDTLL DLQNNKITEI KGDGFKNLKN LHALILVNNK ISKVSPGAF PLVKLERLYL SKNQLKELPE KMPKTLQELR AHENEITKVR KVTFNGLNQM IVIELGTNPL KSSGIENGAF QGMKKLSYIR IADTNITSIP QGLPPSLTEL HLDGNKISRV DAASLKGLNN LAKLGLSFNS ISAVDNGSLA NTPHLRELHL DNNKLTRVPG GLAEHKYIQV VYLHNNNISV VGSSDFCPPG HNTKKASYSG VSLFSNPVQY WEIQPSTFRC VYVRSAILQLG NYK

image



Background

Alternative Names Decorin (Bone proteoglycan II) (PG-S2) (PG40)

Product Name: Recombinant Human DCN Protein
Catalog #: PEH2868



Note

For Research Use Only , Not for Diagnostic Use