


Product Name: Recombinant Human ABCB1 Protein
Catalog #: PEH2873



Summary

Name	MDR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08183
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGHKPDNIKG NLEFRNVHFS YPSRKEVKIL KGLNLKVQSG QTVALVGNSG CGKSTTVQLM QRLYDPTEGM VSVDGQDIRT INVRFLREII GVVSQEPVLF ATTIAENIRY GRENVTMDEI EKAVKEANAY DFIMKLPHKF DTLVGERGAQ LSGGQKQRIA IARALVRNPK ILLDEATSA LDTESEAVVQ VALDKARKGR TTIVIAHRLS TVRNADVIAG FDDGVIVEKG NHDELMKEKG IYFKLVMTQT AGNEVELENA ADESKSEIDA LEMSSNDSRS

image



Background

Alternative Names	ATP-dependent translocase ABCB1 (ATP-binding cassette sub-family B member 1) (Multidrug resistance protein 1) (EC 7.6.2.2) (P-glycoprotein 1) (Phospholipid transporter ABCB1) (EC 7.6.2.1) (CD antigen CD243)
--------------------------	--

Product Name: Recombinant Human ABCB1 Protein
Catalog #: PEH2873



Note

For Research Use Only , Not for Diagnostic Use