


Product Name: Recombinant Human HSP90AB1 Protein
Catalog #: PEH2877



Summary

Name	HS90B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08238
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVETFAFAQE IAQLMSLIIN TFYSNKEIFL RELISNASDA LDKIRYESLT DPSKLD SGKE LKIDIIPNPQ ERTLTLVDTG IGMTKADLIN NLGTIAKSGT KAFMEALQAG ADISMIGQFG VGFYSAYLVA EKVVVITKHN DDEQYAWESS AGGSFTVRAD HGEPIGRGTK VILHLKEDQT EYLEERRVKE VVKKHSQFIG YPITLYLEKE REKEISDDEA EEEKGEKEEE DKDDEEKPKI EDVGSDEEDD SGKDKKKKTK KIKEKYIDQE ELNKT KPIWT RNPDDITQEE YGEFYKSLTN DWEDHLAVKH FSVEGQLEFR ALLFIPRRAP FDLFENK KKK NNIKLYVRRV FIMDSCDELI PEYLN FIRGV VDS EDLPLNI SREMLQQSKI LKVIRKNIVK KCLELFSELA EDKENYKKFY EAFSKNLKLG IHEDSTNRRR LSELLRYHTS QSGDEMTSLS EYVSRMKETQ KSIYYITGES KEQVANS AFV ERVRKRGFEV VYMTEPID EY CVQQ LKEFDG KSLVSVTKEG

image



Background



Product Name: Recombinant Human HSP90AB1 Protein
Catalog #: PEH2877

Alternative Names

Heat shock protein HSP 90-beta (HSP 90) (Heat shock 84 kDa) (HSP 84)
(HSP84)

Note

For Research Use Only , Not for Diagnostic Use