


Product Name: Recombinant Human PTPRC Protein
Catalog #: PEH2881



Summary

Name	PTPRC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08575
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YNKETKLFTA KLVNENVEG GNNTCTNNEV HNLTECKNAS VSISHNSCTA PDKTLILDVP PGVEKFQLHD

image



Background

Alternative Names	Receptor-type tyrosine-protein phosphatase C (EC 3.1.3.48) (Leukocyte common antigen) (L-CA) (T200) (CD antigen CD45)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use