


Product Name: Recombinant Human IL3 Protein
Catalog #: PEH2883



Summary

Name	IL3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08700
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIITHLKQPP LPLDFNNLN GEDQDILMEN NLRPPNLEAF NRAVKSLQNA SAIESILKNL LPCLPLATAA PTRHPIHIKD GDWNEFRRKL

image



Background

Alternative Names Interleukin-3 (IL-3) (Hematopoietic growth factor) (Mast cell growth factor) (MCGF) (Multipotential colony-stimulating factor) (P-cell-stimulating factor)

Note

For Research Use Only , Not for Diagnostic Use