


Product Name: Recombinant Human F7 Protein
Catalog #: PEH2884



Summary

| | |
|------------------------------|---|
| Name | FA7_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P08709 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 18 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | GSLERECKEE QCSFEEAREI FKDAERTKLF WISYSDGDQC ASSPCQNGGS CKDQLQSYIC FCLPAFEGRN CETHKDDQLI CVNENGGCEQ YCSDHTGTKR SCRCHGYSL LADGVSCTPT VEYPCGKIPI LEKRNASKPQ |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Coagulation factor VII (EC 3.4.21.21) (Proconvertin) (Serum prothrombin conversion accelerator) (SPCA) (Eptacog alfa) [Cleaved into: Factor VII light chain; Factor VII heavy chain] |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human F7 Protein
Catalog #: PEH2884

