


Product Name: Recombinant Human KRT19 Protein
Catalog #: PEH2885



Summary

Name	K1C19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08727
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLTMQNLNDR LASYLDKVRA LEAANGELEV KIRDWYQKQG PGPSRDYSHY YTTIQDLRDK ILGATIENSR IVLQIDNARL AADDFR TKFE TEQALRMSVE ADINGLRRVL DELTLARTDL EMQIEGLKEE LAYLKKNHEE EISTLRGQVG GQVSVEVDSA PGTDLAKILS DMRSQYEVMA EQNRKDAEAW FTSRTEELNR EVAGHTEQLQ MSRSEVTDLR RTLQGLEIEL QSQLSMKAAL EDTLAETEAR FGAQLAHIQA LISGIEAQLG DVRADSERQN QEYQRLMDIK SRLEQEIATY RSLLEGQEDH

image



Background

Alternative Names Keratin, type I cytoskeletal 19 (Cytokeratin-19) (CK-19) (Keratin-19) (K19)

Product Name: Recombinant Human KRT19 Protein
Catalog #: PEH2885



Note

For Research Use Only , Not for Diagnostic Use