


Product Name: Recombinant Human MMP7 Protein
Catalog #: PEH2894



Summary

Name	MMP7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P09237
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPLPQEAGGM SELQWEQAQD YLKRFYLYDS ETKNANSLEA KLKEMQKFFG LPITGMLNSR VIEIMQKPRC GVPDVAEYSL FPNPKWTSK VVTYRIVSYT RDLPHITVDR LVSKALNMWG KEIPLHFRKV VWGTADIMIG FARGAHGDSY PFDGPGNTLA HAFAPGTGLG GDAHFDEDER WTDGSSLGIN FLYAATHELG HSLGMGHSSD PNAVMYPTYG NGDPQNFKLS QDDIKGIQKL YGKRSNSRKK

image



Background

Alternative Names	Matrilysin (EC 3.4.24.23) (Matrin) (Matrix metalloproteinase-7) (MMP-7) (Pump-1 protease) (Uterine metalloproteinase)
--------------------------	---

Note

Product Name: Recombinant Human MMP7 Protein
Catalog #: PEH2894



For Research Use Only , Not for Diagnostic Use