


Product Name: Recombinant Human HMGB1 Protein
Catalog #: PEH2897



Summary

Name	HMGB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09429
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGKGDPPKPR GKMSYAFFV QTCREEHKKK HPDASVNFSE FSKKCSERWK TMSAKEKGKF EDMAKADKAR YEREMKTYIP PKGETKKKFK DPNAPKRPPS AFFLFCSEYR PKIKGEHPGL SIGDVAKKLG EMWNNTAADD KQPYEKKA LKEKYEKIDIA AYRAKGKPD AAKGVVKA EK SKKKKEEEED EEEEEDEEE EDEEDEDEEE DDDDE

image



Background

Alternative Names High mobility group protein B1 (High mobility group protein 1) (HMG-1)

Note

Product Name: Recombinant Human HMGB1 Protein
Catalog #: PEH2897



For Research Use Only , Not for Diagnostic Use