


Product Name: Recombinant Human PARP1 Protein
Catalog #: PEH2900



Summary

Name	PARP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P09874
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GHSIRHPDVE VDGFSLELRWD DQQKVKKTAE AGGVTGKGQD GIGSKAEKTL GDFAAEYAKS NRSTCKGCME KIEKGQVRLS KKMVDPEKPQ LGMIDRWYHP GCFVKNREEL GFRPEYSASQ LKGFSLATE DKEALKKQLP GVKSEGKRKG DEV DGVDEVA KKKSKKEKDK DSKLEKALKA QNDLIWNIKD ELKKVCSTND LKELLIFNKQ QVPSGESAIL DRVADGMVFG ALLPCEEC SG QLVFKSDAYY CTGDVTAWTK CMVKTQTPNR KEWVTPKEFR EISYLKKLV KKQDRIFPPE TSASVAATPP

image



Background

Alternative Names	Poly [ADP-ribose] polymerase 1 (PARP-1) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 1) (ARTD1) (DNA ADP-
--------------------------	--

Product Name: Recombinant Human PARP1 Protein
Catalog #: PEH2900



ribosyltransferase PARP1) (EC 2.4.2.-) (NAD(+) ADP-ribosyltransferase 1)
(ADPRT 1) (Poly[ADP-ribose] synthase 1) (Protein poly-ADP-
ribosyltransferase PARP1) (EC 2.4.2.-)

Note

For Research Use Only , Not for Diagnostic Use