


Product Name: Recombinant Human CSF3 Protein
Catalog #: PEH2901



Summary

| | |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name | CSF3_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-MBP & C-His |
| Accession # | P09919 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 61 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | QSFLKCLEQ VRKIQGDGAA LQEKLVSECA TYKLCHPEEL VLLGHSLGIP WAPLSSCPSQ ALQLAGCLSQ LHSGLFLYQG LLQALEGISP ELGPTLDTLQ LDVADFATTI WQQMEELGMA PALQPTQGAM PAFASAFQRR AGGVLVASHL QSFLEVSYRV |

image



Background

Alternative Names Granulocyte colony-stimulating factor (G-CSF) (Pluripoietin) (Filgrastim) (Lenograstim)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CSF3 Protein
Catalog #: PEH2901

