



**Product Name: Recombinant Human CHGA Protein**  
**Catalog #: PEH2913**



## Summary

<b>Name</b>	CMGA_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P10645
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	17 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RLSKEWEDSK RWSKMDQLAK ELTAEKRLEG QEEEEENRDS SMKLSFRARA YGFRGPGPQL RRGWRPSSRE DSLEAGLPLQ VRGYPEEKKE EEGSANRRPE DQELESLSAI EAELEKVAHQ

## image



## Background

**Alternative Names** Chromogranin-A (CgA) (Pituitary secretory protein I) (SP-I) [Cleaved into: Vasostatin-1 (Vasostatin I); Vasostatin-2 (Vasostatin II); EA-92; ES-43; Pancreastatin; SS-18; WA-8; WE-14; LF-19; Catestatin (SL21); AL-11; GV-19; GR-44; ER-37; GE-25; Serpinin-RRG; Serpinin; p-Glu serpinin precursor]

## Note

**Product Name: Recombinant Human CHGA Protein**  
**Catalog #: PEH2913**



---

For Research Use Only , Not for Diagnostic Use