


Summary

Name	ITAM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P11215
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQQPQKFPEA LRGCPQEDSD IAFLIDGSGS IIPHDFRRMK EfvSTVMEQL KKSktLfSLM QYSEEFRIHF TFKEFQNNPN PRSLVKPITQ LLGRtHTATG IRKVVRELFN ITNGARKNAF KILVVITDGE KFGDPLGYED VIPEADREGV IRYVIGVGDA FRSEKSRQEL NTIASKPPRD HVFQVNNFEA LKTIQNQLRE KIFAIEGTQT GSSSSFEHEM

image



Background

Alternative Names	Integrin alpha-M (CD11 antigen-like family member B) (CR-3 alpha chain) (Cell surface glycoprotein MAC-1 subunit alpha) (Leukocyte adhesion receptor MO1) (Neutrophil adherence receptor) (CD antigen CD11b)
--------------------------	--

Product Name: Recombinant Human ITGAM Protein
Catalog #: PEH2918



Note

For Research Use Only , Not for Diagnostic Use