


**Product Name: Recombinant Human CYP2C9 Protein**  
**Catalog #: PEH2922**



## Summary

<b>Name</b>	CP2C9_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P11712
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	10 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	HQESMDMNNP QDFIDCFLMK MEKEKHNQPS EFTIESLENT AVDLFGAGTE TTSTTLRYAL LLLLKHPEVT

## image



## Background

**Alternative Names** Cytochrome P450 2C9 (EC 1.14.14.1) ((R)-limonene 6-monooxygenase) (EC 1.14.14.53) ((S)-limonene 6-monooxygenase) (EC 1.14.14.51) ((S)-limonene 7-monooxygenase) (EC 1.14.14.52) (CYPIIC9) (Cholesterol 25-hydroxylase) (Cytochrome P-450MP) (Cytochrome P450 MP-4) (Cytochrome P450 MP-8) (Cytochrome P450 PB-1) (S-mephenytoin 4-hydroxylase)

## Note

**Product Name: Recombinant Human CYP2C9 Protein**  
**Catalog #: PEH2922**



---

For Research Use Only , Not for Diagnostic Use