


Product Name: Recombinant Human KRT3 Protein
Catalog #: PEH2924



Summary

Name	K2C3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12035
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QIKTLNNKFA SFIDKVRFLE QQNKVLETKW NLLQQQGTSS ISGTNNLEPL FENHINYLRSL YLDNILGERG RLDSELKNME DLVEDFKKKY EDEINKRTAA ENEFVTLKKD VDSAYMKNKE LQAKVDALID EIDFLRTLTD AELSQQMSHI SDTSVVLSDM NNRLDLDSI IAEVRAQYED IAQRSKAEAE ALYQTKLGEL QTTAGRHGDD LRNTKSEIIE LNRMIQRLRA EIEGVKKQNA NLQTAIAEAE QHGEMLKDA NAKLQELQAA LQQAQDDLAR LLRDYQELMN VKLALDVEIA TYRKLLEGEE YRMSGECPSA VSISVVSST

image



Background

Alternative Names	Keratin, type II cytoskeletal 3 (65 kDa cytokeratin) (Cytokeratin-3) (CK-3) (Keratin-3) (K3) (Type-II keratin Kb3)
--------------------------	--

Product Name: Recombinant Human KRT3 Protein
Catalog #: PEH2924



Note

For Research Use Only , Not for Diagnostic Use