


Product Name: Recombinant Human HARS1 Protein
Catalog #: PEH2925



Summary

Name	HARS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12081
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FR TICSSVDK LDKVSWEEVK NEMVGEKGLA PEVADRIGDY VQQHGGVSLV EQLLQDPKLS QNKQALEGLG DLKLLFEYLT LFGIDDKISF

image



Background

Alternative Names	Histidine--tRNA ligase, cytoplasmic (EC 6.1.1.21) (Histidyl-tRNA synthetase) (HisRS)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use