

## Summary

<b>Name</b>	PLSL_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P13796
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	62 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RKAINKKEGI CAIGGTSEQS SVGTQHSYSE EEKYAFVNW I NKALENDPDC RHVIPMNPNT NDLFNAVGDG IVLCKMINLS VPDTIDERTI NKKKLTPTFI QENLNLALNS ASAIGCHVVN IGAEDLKEGK PYLVLGLLWQ VIKIGLFADI ELSRNEALIA LLREGESLED LMKLSPEELL LRWANYHLEN AGCNKIGNFS TDIKDSKAYY HLLEQVAPKG DEEGVPAVVI DMSGLREKDD IQRAECMLQQ AERLGCRQFV TATDVVRGNP KLNLAFIANL FNRYPALHKP ENQDIDWGAL EGETREERTF RNWMNSLGVN PRVNHLYSDL SDALVIFQLY EKIKVPVDWN RVNKPPYPKL GGNMKKLENC NYAVELGKNQ AKFSLVGIGG QDLNEGNRTL TLALIWQLMR RYTLNILEEI GGGQKVNDI IVNWVNETLR EAKKSSSISS FKDPKISTSL PVLDLIDAIQ PGSINYDLLK TENLNDDEKL NNAKYAISMA RKIGARVYAL PEDLVEVNP K MVMTVFACLM

## image



## Background

**Product Name: Recombinant Human LCP1 Protein**  
**Catalog #: PEH2930**



---

**Alternative Names**

Plastin-2 (L-plastin) (LC64P) (Lymphocyte cytosolic protein 1) (LCP-1)

**Note**

For Research Use Only , Not for Diagnostic Use