


Product Name: Recombinant Human HSP90B1 Protein
Catalog #: PEH2931



Summary

Name	ENPL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P14625
Host	E.coli
Species	Human
Molecular Mass	68 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	REEEAIQLDG LNASQIRELR EKSEKFAFQA EVNRMMKLII NSLYKNKEIF LRELISNASD ALDKIRLISL TDENALSGNE ELTVKIKCDK EKNLLHVTD GVGMTREELV KNLGTIAKSG TSEFLNKMTE AQEDGQSTSE LIGQFGVGFY SAFLVADKVI VTSKHNDTQ HIWESDSNEF SVIADPRGNT LGRGTTITLV LKEEASDYLE LDTIKNLVKK YSQFINFPIY VWSSKTETVE EPMEEEEAAK EEKEESDDEA AVEEEEEEEK PKTKKVEKTV WDWELMNDIK PIWQRPSKEV EEDEYKAFYK SFSKESDDPM AYIHFTAEGE VTFKSILFVP TSAPRGLFDE YGSKKS DYIK LYVRRVITD DFHDMMPKYL NFVKGVVDS DLPLNVSRET LQQHKLLKVI RKKLVRKTLD MIKKIADDKY NDTFWKEFGT NIKLGVIEDH SNRTRLAKLL RFQSSHPTD ITSLDQYVER MKEKQDKIYF MAGSSRKEAE SSPFVERLLK KGYEVIYLTE PVDEYCIQAL PEFDGKRFQN VAKEGVKFDE SEKTKESREA

image



Product Name: Recombinant Human HSP90B1 Protein
Catalog #: PEH2931



Background

Alternative Names

Endoplasmin (94 kDa glucose-regulated protein) (GRP-94) (Heat shock protein 90 kDa beta member 1) (Tumor rejection antigen 1) (gp96 homolog)

Note

For Research Use Only , Not for Diagnostic Use