


Product Name: Recombinant Human CCNB1 Protein
Catalog #: PEH2932



Summary

Name	CCNB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P14635
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEPEPEPEPV KEEKLSPEPI LVDTASPSM ETSGCAPAEE DLCQAFSDVI LAVNDVDAED GADPNLCSEY VKDIYAYLRQ LEEEQAVRPK YLLGREVTGN MRAILIDWLV QVQMKFRLQ

image



Background

Alternative Names G2/mitotic-specific cyclin-B1

Note

For Research Use Only , Not for Diagnostic Use