



Product Name: Recombinant Human EZR Protein
Catalog #: PEH2940



Summary

Name	EZRI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15311
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVQQMKAQAR EEKHQQLER QQLETEKKRR ETVEREKEQM MREKEELMLR LQDYEEKTKK AERELSEQIQ RALQLEEERK RAQEEAERLE ADRMAALRAK EELERQAVDQ IKSQEQLAAE LAEYAKIAL LEEARRRKED EVEEWQHRAK EAQDDLKTK EELHLVMTAP PVPVYEP VSYHVQESLQ DEGAEPTGYS AELSSEGIRD DRNEEKRITE AEKNERVQRQ LLLSSELSQ ARDENKRTHN DIIHNENMRQ GRDKYKTLRQ IRQGNTKQRI

image



Background

Alternative Names Ezrin (Cytovillin) (Villin-2) (p81)

Product Name: Recombinant Human EZR Protein
Catalog #: PEH2940



Note

For Research Use Only , Not for Diagnostic Use