


Product Name: Recombinant Human ELN Protein
Catalog #: PEH2942



Summary

Name	ELN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15502
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPGVGLPGVY PGGVLPGARF PGVGVLPGVP TGAGVKPKAP GVGGAFAGIP GVGPFGGPQP GVPLGYPIKA PKLPGGYGLP YTTGKLPGY GPGGVAGAAG KAGYPTGTGV GPQAAAAAAAA KAAAKFGAGA AGVLPGVGGA GVPGVPGAIP GIGGIAGVGT PAAAAAAAAA AKAAYGAAA GLVPGGPGFG PGVGVPGAG VPGVGVPGAG IPVVPAGIP GAAVPGVVSP EAAAKAAKA AKYGARPGVG VGGIPTYGVG AGGFPGFGVG VGGIPGVAGV PGVGGVPGVG GVPGVGISPE AQAAAAAKAA KYGAAGAGVL GGLVPGAPGA VPGVPGTGGV PGVGTAAAA AKAAAKAAQF GLVPGVGVAP GVGVPAGVG APGVGLAPGV GVAPGVGVAP

image



Background

Alternative Names	Elastin (Tropoelastin)
--------------------------	------------------------

Product Name: Recombinant Human ELN Protein
Catalog #: PEH2942



Note

For Research Use Only , Not for Diagnostic Use