


Product Name: Recombinant Human PDGFRA Protein
Catalog #: PEH2947



Summary

| | |
|------------------------------|--|
| Name | PGFRA_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P16234 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 11 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | SYEKIHLDFL KSDHPAVARM RVDSDNAYIG VTYKNEEDKL KDWEGGLDEQ RLSADSGYII PLPDIDPVPE EEDLGKRRH |

image



Background

| | |
|--------------------------|---|
| Alternative Names | Platelet-derived growth factor receptor alpha (PDGF-R-alpha) (PDGFR-alpha) (EC 2.7.10.1) (Alpha platelet-derived growth factor receptor) (Alpha-type platelet-derived growth factor receptor) (CD140 antigen-like family member A) (CD140a antigen) (Platelet-derived growth factor alpha receptor) (Platelet-derived growth factor receptor 2) (PDGFR-2) (CD antigen CD140a) |
|--------------------------|---|

Note

Product Name: Recombinant Human PDGFRA Protein
Catalog #: PEH2947



For Research Use Only , Not for Diagnostic Use