


Product Name: Recombinant Human POR Protein
Catalog #: PEH2951



Summary

Name	NCPR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P16435
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSSVRESSFV EKMKKTGRNI IVFYGSQTGT AEEFANRLSK DAHRYGMRGM SADPEEYDLA DLSSLPEIDN ALVVFCMATY GEGDPTDNAQ DFYDWLQETD VDLSGVKFAV FGLGNKTYEH FNAMGKYVDK RLEQLGAQRI FELGLGDDDG NLEEDFITWR EQFWPAVCEH FGVEATGEES SIRQYELVH

image



Background

Alternative Names NADPH--cytochrome P450 reductase (CPR) (P450R) (EC 1.6.2.4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human POR Protein
Catalog #: PEH2951

