


Product Name: Recombinant Human IGFBP2 protein
Catalog #: PEH2961



Summary

Name	IBP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P18065
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CEKRRDAEYG ASPEQVADNG DDHSEGGGLVE NHVDSTMNML GGGGSAGRKP LKSGMKELAV FREKVTEQHR QMGKGGKHHL GLEPKKLRP PPARTPCQQE LDQVLERIST MRLPDERGPL EHLVSLHIPN CDKHGLYNLK QCKMSLNGQR GECWCVPNPNT GKLIQGAPTI RGDPECHLFY NEQQEARGVH

image



Background

Alternative Names	Insulin-like growth factor-binding protein 2 (IBP-2) (IGF-binding protein 2) (IGFBP-2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IGFBP2 protein
Catalog #: PEH2961

