


Product Name: Recombinant Human EIF2S2 Protein
Catalog #: PEH2971



Summary

Name	IF2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20042
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IDEAEEGVKD LKIESDVQEP TEPEDDLDIM LGNKKKKKKKN VKFPDEDEIL EKDEALEDED NKKDDGISFS NQTGPAWAGS ERDYTYEELL NRVFNIMREK NPDMVAGEKR KFVMKPPQVV RVGTTKTSFV NFTDICKLLH

image



Background

Alternative Names	Eukaryotic translation initiation factor 2 subunit 2 (Eukaryotic translation initiation factor 2 subunit beta) (eIF-2-beta)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use